

Temporary pacing leads



Reliable pacing. Precise placement. Easy handling.

A range of temporary pacing leads and balloon-tipped catheters for dependable performance when it matters most.

Available
through
NHSSC

Key features

- Non-braided, atraumatic shaft for smooth insertion
- Pre-shaped curves for optimal positioning
- Stainless steel electrodes for reliable signal
- 4, 5, 6 French diameters
- Sterile, disposable, latex-free
- Shrouded and non-shrouded options
- Additional tip styles available on request



SPIKE LC bipolar temporary pacing leads

Designed for smooth, controlled positioning, SPIKE LC pacing leads provide dependable pacing and recording performance. Available in 4F, 5F and 6F sizes with shrouded or non-shrouded pins.



SPIKE FLOW balloon-tipped pacing catheter

The SPIKE FLOW bipolar balloon catheter offers flow-directed placement and can be positioned without fluoroscopy when required. Features depth markers at 10cm intervals to aid placement and procedural efficiency.

Bipolar Temporary Pacing Leads

	Code	Size	Electrode spacing	Usable length
Spike LC Pacing Lead with shrouded pins	52204S	4F	10mm	90cm
	52205S	5F	10mm	125cm
	52206S	6F	10mm	125cm
Spike LC Pacing Lead with non-shrouded pins, with protective caps	52204N	4F	10mm	90cm
	52205N	5F	10mm	125cm
	52206N	6F	10mm	125cm
Spike LC Pacing Lead, balloon tipped with non-shrouded with protective caps	12165	5F	10mm	110cm

SPIKE™ is a registered trademark of FIAB

All trademarks are acknowledged as the property of their respective owners.

Get in touch:

+44 (0)1606 331119 | info@apccv.com | apccardiovascular.co.uk/contact

TPL1026

Scan for more
information

